

Introduction

Nondestructive Infra Red Light High-throughput Biochemical Analysis of Single Seed Canola

Daniel Gachotte
Randy Pell
Paul Morabito
Ted Weglarz

1/14/2008

Introduction

Omega-9 Canola oil is the ideal
product for the rapidly growing
new segment of healthy,
naturally-stable vegetable oils



Trademark of Dow AgroSciences LLC

Introduction


Omega-9 Canola Oil Products Match Consumer Need

- Food industry seeking oxidative stability have used partially hydrogenated vegetable oils
 - Eight billion pounds of partially hydrogenated oils used in U.S. (overall market value approximately \$2.8B)
- FDA announces mandatory labeling of trans fats (Jan. 06')
- Omega-9* canola increases content of “good” fats and reduces “bad” fats, while providing critical functionality
- Decreased risk of coronary heart disease associated with “bad fats”

1/14/2008

*Trademark of Dow AgroSciences LLC

Growth of Dow AgroSciences Omega-9 Oils Program



1/14/2008

Introduction

Summary

- **Automation: hardware and software**
- **Chemometrics: model prediction**
- **Improvements: new instrumentation**
- **Application: large population screening**

1/14/2008

Automation

Generation of canola single seed NIR spectra



OPUS software driving data acquisition

1/14/2008

Automation

Generation of canola single seed NIR spectra



1/14/2008

Automation

Efficiency achieved in single seed


- **Per Tecan RSP150 in 24 hours period**
 - **throughput increased by 35%** (1800 to 2250 seed from envelopes)
 - **throughput increased by 50%** (2880 seeds from plates)
 - **missing seed decreased using new vacuum chip** (1% from plate to plate to 2% from envelope)

1/14/2008

Automation

Software/Hardware Interactions

- Gemini controls robot
- VB program triggered by Gemini
- VB program triggers Opus
- Opus scans spectrum
- InStep™ triggered by VB
- Gemini resumes when VB complete



1/14/2008

Summary

- **Automation: hardware and software**
- **Chemometrics: model prediction**
- **Improvements: new instrumentation**
- **Application: large population screening**

1/14/2008

Chemometrics


Calibration Model Development

- Plot the data
 - Identify obvious problem data
- Construct model
 - Select a preprocessing method
 - Select samples (outlier detection)
 - Select variables (wavelength selection)
- Validate model
 - Examine calibration diagnostics
 - Assess performance
 - cross validation often used
 - separate test set prediction

1/14/2008

Chemometrics

Select preprocessing




1/14/2008

Chemometrics

2nd derivative of NIR data

25 point window, 3rd order polynomial



1/14/2008

Chemometrics

Validation

Cross validation

	Sample # 1	Sample # 2	Sample # 3	Sample # 4	Sample # 5	Sample # 6	Sample # 7	Sample # 8
Cycle #1	v	c	c	c	c	c	c	c
Cycle #2	c	v	c	c	c	c	c	c
Cycle #3	c	c	v	c	c	c	c	c
Cycle #4	c	c	c	v	c	c	c	c
Cycle #5	c	c	c	c	v	c	c	c
Cycle #6	c	c	c	c	c	v	c	c
Cycle #7	c	c	c	c	c	c	v	c
Cycle #8	c	c	c	c	c	c	c	v

- This is a graphic example of a LOO-CV experiment.
- Eight samples lead to 8 calibration and prediction cycles
- Each sample left out (put in validation) once and only once.
- Should be possible to execute in any reasonable software package.

1/14/2008

Summary

- **Automation: hardware and software**
- **Chemometrics: model prediction**
- **Improvements: new instrumentation**
- **Application: large population screening**

1/14/2008

Improvements

Acquisition of new FTNIR: MPA

- Two fiber optic probe
- Multi purpose analyzer (transmission and reflection)
- New version of OPUS
- Same rugged interferometer




1/14/2008

Improvements

System modification

- New integrated sphere (8% intensity gain)
- High intensity source (2.5 X signal amplitude)



- Dry N₂ purged instrument to eliminate residual moisture increase reliability and stability
- New version of OPUS (firmware update)

1/14/2008

Summary

- **Automation: hardware and software**
- **Chemometrics: model prediction**
- **Improvements: new instrumentation**
- **Application: large population screening**

1/14/2008

Applications

Generation of canola single seed reference methods


- Quantitative and qualitative analysis require
 - accurate weighing (0.01 mg)
 - exhaustive extraction (multiple)
 - automated pipetting (Hamilton)
- Analytical instrument
 - GC/FID
 - HPLC/Fluorescence/ELSD/UV
- Sensitivity and reproducibility
 - Medium throughput (96 well plates)
 - Derivatization to increase sensitivity

1/14/2008

Applications

Prediction Performance for oil

Factor	RMSEP	SEP	r2	Range/RMSEP	SD/RMSEP
4*	0.0483	0.0483	0.9679	26.3	5.61
3*	1.0429	1.0429	0.9138	14.02	3.42
6*	1.53	1.5299	0.9588	21.09	4.95





1/14/2008

Applications

Prediction Performance for Protein


Factor	RMSEP	SEP	r2	Range/RMSEP	SD/RMSEP
4*	0.0275	0.028	0.966	21.6	5.49
3*	0.7836	0.783	0.884	12.62	2.96




1/14/2008

Applications

Single seed screening process



1/14/2008



Acknowledgements

TDCC

Jon Zieman
Rick Winterton
Bryant LaFreniere

External companies

Tecan GMBH
QTA (Cognis)
Controldevelopment
Brucker Optics

DAS

Steve Thompson
Chuck Kahl
Tom Patterson

Max Ruegger
Sam Reddy
Larry Sernyk

Jana Alderson
Beth Blakeslee
Scott Greenwalt

Thank you for your attention